Canadian Light Source Soft X-ray Spectromicroscopy (SM) Beamline Manual Scanning Transmission X-ray Microscopy (STXM) Branch



©

Canadian Light Source Inc. University of Saskatchewan 101 Perimeter Road Saskatoon, Saskatchewan Canada S7N 0X4 SM beamline Phone: 1-306-657-3609 CLS Main Phone: 1-306-657-3500 CLS User Services Phone: 1-306-657-3700 Fax: 1-306-657-3535 CLS website: http://www.lightsource.ca/

10ID-1 CLS SM Beamline Manual

Version: 13 June 2011

Table of Contents

Table of Acronyms	5
Scope of the Document	7
SM Beamline Staff and Beam Team Leaders	9
1.0 Check List before Your Arrival	. 11
1.1 Confirm SM beamline time	. 11
1.2 Perform Proposal Amendment if Required	. 11
1.3 CLS Safety Training	. 11
1.4 Operation Manuals and Quality Assurance Documentation	. 11
2.0 At the SM Beamline	. 13
2.1 SM Beamline Specific Orientation (BSO)	. 13
2.2 Permit & Pink Sheet	. 13
2.3 Food and Drink	. 13
2.4 Floor Coordinator	. 13
2.5 Invoice	. 13
2.6 CLS Beamline Staff Availability	. 15
3.0 Overview of the SM Beamline and STXM Microscope	. 17
4.0 PART A: Enabling the SM Beamline	. 19
4.1 Starting the STXM Branch Graphical User Interface	. 19
4.2 Branch Selection	. 19
4.3 Safety and Photon Shutters/Access Control Information System Panel	. 19
4.4 Beamline Safety and SM-Photon Shutters and the Variable Aperture	. 23
4.5 Closing the SM Beamline (at the end of the shift)	. 23
4.6 Enabling the SM Beamline (Quick Guide)	. 25
5.0 PART B: STXM Microscope Set-up	. 27
5.1 Zone Plate Selection	. 27
5.2 Grating Selection	. 27
5.3 Energy Selection	. 31
5.4 Harmonic Selection	. 33
5.5 Polarization Selection	. 33
5.6 Constant of Fixed Focus (Cff)	. 35
6.0 PART C: STXM Microscope Optimization	. 37
6.1 Verify that Beam is reaching the Detector	. 37
6.2 Detector Scan	. 39
6.3 OSA Scan	. 41
6.4 OSA Focusing Scan	. 41
6.5 Exit Slits	. 43
6.6 M3STXM Mirror Pitch	. 43
6.7 EPU Offset	. 45
6.8 STXM Microscope Set-up and Optimization Summary	. 47
7.0 SM Beamline (STXM) Energy Calibration	. 49
Appendix A: Selecting Parameters for Specific Elements	. 51
Appendix B: LOADING/UNLOADING a Sample into the STXM TANK	. 53
Appendix C: SM Beamline Polarization Setup	. 55
Appendix D: Zone Plate, Order Sorting Aperture (OSA), Sample & Detector Relationship	. 57
Appendix E: SM Data Access and Data Transfer	. 59

Appendix F: Connecting to the Printer in the STXM hutch	63
Appendix G: SM Beamline – Computer's and Software for Users	65
Appendix H: Selecting the optimum Cff value	67
Appendix I: References	69

Filename: 91013_SM_BeamlineManual.doc

Acronym	Meaning			
CLS	Canadian Light Source			
SM	Spectromicroscopy			
STXM	Scanning Transmission X-ray Microscope			
PEEM	X-ray Photoemission Electron Microscope			
ACIS	Access Control Information System			
SM_BEAMLINE_GUI	User control panel for SM Beamline			
STXM_Control	Name of STXM control program			
STXM_CONTROL_GUI	User control panel for STXM			
Lox	Name of PEEM control program			
cff	Constant of fixed focus			
α/β	mirror and grating angles			

Table of Acronyms

Scope of the Document

The purpose of this manual is to help novice and less experienced Users to set-up the SM beamline and set-up and optimize the STXM End station (i.e., the microscope), in a step by step procedure, including snap shots of the actual buttons to push. Many of the steps are discussed in detail, giving the User some insight into what they are doing, and what the expected outcome should be.

After following the steps in this manual the User will be able to enable the SM beamline (Part A) and set-up (Part B) and optimized (Part C) the STXM Microscope, setting the stage for Data Acquisition. For detailed discussion on Data Acquisition see the STXM USER MANUAL. Note that the SM Beamline needs to be enabled first prior to setting up the STXM microscope. Also, is included a section on gas calibration.

There is also an Appendix section that deals with Sample Load/Unloading, and Zone plate and sample relationship.

This manual also provides Users with some useful information, including the SM Beamline staff contact information, what a User needs to check before their arrival and what the User can expect upon their arrival at the SM Beamline.

SM Beamline Staff and Beam Team Leaders



Chithra Karunakaran Staff Scientist Tel: 1-306-657-3749 chithra.karunakaran@lightsource.ca



Jian Wang Staff Scientist Tel: 1-306-657-3546 jian.wang@lightsource.ca



Uday Lanke PEEM Research Associate Tel: 1-306-657-3769 uday.lanke@lightsource.ca



Yingshen Lu Science Associate Tel: 1-306-657-3743 Yingshen.lu@lightsource.ca



James Dynes STXM Research Associate Tel: 1-306-657-3840 james.dynes@lightsource.ca



Adam Hitchcock Beam Team Leader, McMaster University Tel: 1-905-525-9140 Ext. 24749 <u>aph@mcmaster.ca</u> web: http://unicorn.mcmaster.ca



Stephen Urquhart ity PEEM Leader, University of Saskatchewan Tel: 1-306-966-4657 <u>stephen.urquhart@usask.ca</u> web: http://homepage.usask.ca/~sgu703/index.html

10ID-1 CLS SM Beamline Manual

Version: 13 June 2011

1.0 Check List before Your Arrival

1.1 Confirm SM beamline time

The CLS User Services office usually sends out a remainder to Users of their upcoming beam time one week prior to beam time. The most current beamline schedule is available on the CLS website (http://www.lightsource.ca/uso/) or ask the beamline staff for a copy. If you have any questions or concerns, please do not hesitate to contact the beamline staff or the CLS User Services office.

1.2 Perform Proposal Amendment if Required

It is the responsibility of the User to make sure that proposal is up to date, particularly that the samples they intend to run during their shifts at the SM beamline are all listed in their proposal. If the User plans to run samples not currently listed in their proposal, the proposal needs to be amended prior to their beam time and Health, Safety & Environment (HSE) must approve the changes. The amendments must be done on-line, editing the previously submitted proposal. In order to ensure that HSE gives their approval, the changes should be made at least 4 weeks prior to their beam time. Samples not listed in the proposal will not be allowed to be run at the SM beamline until the SM Beamline Staff receives HSE approval. Contact the CLS User Services office if you have any questions regarding proposal amendment.

1.3 CLS Safety Training

All Users must have the appropriate safety training prior to operating any beamline or using the laboratories. The training modules can be reviewed and the exams for each module completed on-line or at the time of registration. Arrangements to complete the training at the CLS should be made with the User Services office prior to your arrival. The modules include the Health & Safety Orientation (HSO), Radiation Awareness Module (RAM), Workplace Hazardous Materials Information System (WHMIS) and Laboratory Safety. If radiological materials are to be used, the User must complete the General Radiological Training (GRT). All training expires after 2 years, except WHMIS which does not expire. The CLS User Services office hours are Monday to Friday from 8:00 a.m. to 4:30 p.m.

1.4 Operation Manuals and Quality Assurance Documentation

The SM beamline and microscope STXM operational manuals are available at the SM beamline website (http://exshare.lightsource.ca/sm/). Read them before your arrival at the CLS/beam time so that you can efficiently use your beam time while at the SM beamline. The SM beamline also has a quality assurance document which states the expected performance of the beamline under set conditions.

Spectromicroscopy Beamline Quality Assurance Checklist

Zone Plate			
OSA Hole Size (um) 40 50	60 70	80	
Detector Scan Completed			
OSA Scan Completed			
OSA Focus Scan Completed _			
Exit Slits			
Counts at 395 eV	Hz/s		
Other Comments			
Beamline Staff that set-up and op	timized the ST2	XM Microscope	

Print Name

Signature

Date

Figure 2.1. Spectromicroscopy Beamline Quality Assurance Checklist.

BEAMLINE UNATTENDED	BEAMLINE NOT IN USE
Person to Contact in Case of EMERGENCY:	Date & time you left the beamline:
Contact phone number and/or location:	Notes:
Date & time you left the beamline:	
Expected date & time at which you will return to the beamline:	
Would you like to be contacted if the beam goes down? Yes No	
Notes:	
This form must be posted at the Health & Safety Information Centre of the beamline if you plan to be away more than <u>30 minutes</u> while an experiment is running.	This form must be posted at the Health & Safety Information Centre of the beamline if you stop running your experiment before the end of your scheduled shift. Please notify the Floor Coordinator at 657-3639 prior to leaving.

Figure 2.2. Information sheet (pink sheet) to be left with the permit when the User is away from the beamline for more than 30 mins.

2.0 At the SM Beamline

2.1 SM Beamline Specific Orientation (BSO)

After successfully completing the CLS safety training the SM beamline staff will provide the User with instructions on the procedures and practices Users (i.e. Beamline Specific Orientation (BSO)) must follow to safely operate and protect the User, beamline and microscope. The beamline staff will complete the BSO checklist, and both the staff and User must sign the BSO checklist before the User will be permitted to operate the beamline. The BSO must be completed every 2 years. Also, the SM Beamline Staff will provide the User with a Quality assurance checklist that shows that the STXM microscope is performing as expected (Figure 2.1).

The Beamline Staff will ensure that the microscope is set-up initially for the Users first shift. For first time Users or those needing a refresher, the Beamline Staff will go through all the steps for initially setting up and optimizing the microscope during the first shift and any subsequent shifts as needed.

2.2 Permit & Pink Sheet

Prior to using the beamline, the permit must be signed by both a beamline Staff and a User listed on the permit. Make sure that all the samples that are going to be examined are listed and that all approvals have been signed off (e.g., BSO).

If you leave the beamline for more than 30 minutes you must fill out the HSE form **BEAMLINE UNATTENDED/BEAMLINE NOT IN USE** (i.e., pink sheet) (Figure 2.2) and place it with your permit. The purpose of the pink sheet is so that the CLS Emergency Response Team is made aware of who is on-site in case of an Emergency (e.g., fire).

2.3 Food and Drink

Food and drink are allowed at the SM beamline in the User work area. Only drinks can be taken into the STXM hutch and to the PEEM end station.

2.4 Floor Coordinator

The Floor Coordinator's (FC) is the primary contact for Users while at the SM Beamline, acting as liaison between the Controls Operator, Users and HSE personnel. Their roles are to communicate CLS Policies/Procedures to Users, monitor compliance with those policies, including inspection of the Experiment Permit and perform periodic inspections of all experimental areas. They provide first response in the event of accidents/incidents on the experimental floor. They can enable/disable all beamlines. For FC assistance call 657-**3639**.

2.5 Invoice

In your folder will be an invoice for your shifts. We encourage you to pay while at the CLS. Remit payment at the CLS Reception (1-306-657-3500) during normal office hours (8am to 4:30pm) or to the Floor Coordinator outside of normal office hours. Payment methods include: Cash, Cheque, USASK CFOAPAL, Credit Card (Visa, Mastercard). Make cheques payable to "Canadian Light Source Inc." If you have any questions please contact the CLS User Services at 1-306-657-3700 or <u>clsuo@lightsource.ca</u>.

2.6 CLS Beamline Staff Availability

This page was taken from document 6.7.1.1 Rev. B, dated 2009-02-03. Please see the complete document for more information.

Role of the beamline staff:

Beamline staff will set up the beamline. They will be present at the scheduled start of the experiment to complete the necessary paperwork, and provide any required safety training. They will train the users to operate the beamline, and will assist new user groups through the first few hours of their experiment.

The beamline staff will assist users with their experiments. However, it is not their role to perform experiments for users, including Beam Team members, except by special arrangement. It is the users' responsibility to provide enough qualified individuals to run the experiment for the duration of the beamtime – note that the CLS has 24 hour operation.

Normal working hours:

During normal working hours (9am-5pm Mon-Fri) the beamline staff will be present at the CLS, or able to respond at short notice. Note that under the present schedule only three or four user shifts per week fall within these times.

Outside normal working hours:

Up to 10pm Mon-Fri, and 10am to 10pm weekends and holidays, users may contact the beamline staff for assistance using the contact details provided. If telephoning is indicated, users should use the beamline phone, so that staff can identify that this is likely to be a request for user assistance. A response is not guaranteed, but beamline staff will endeavour to inform users, in advance, of times when no coverage is available. Staff may attempt to solve problems remotely, but may come in to the lab, at their discretion.

Outside these hours we request that users do not attempt to contact the beamline staff, unless prior arrangements have been made for additional support.

In the event that the Floor Co-ordinator (FC) is familiar with the beamline, they may be able to offer assistance. However, in general the FC is not authorised to perform technical fixes, and users should not expect the FC to solve technical problems.

Since 24 hour support cannot be provided, users should be prepared that on occasion no help will be available until the start of normal working hours.

Additional support:

Support beyond that detailed above may be available by prior agreement between the users and the beamline staff. As much as possible, arrangements will be make to provide the highest level of assistance for new user groups. If significant assistance with experimental planning, data acquisition and interpretation is required, it is normally expected that the users will treat the beamline staff as collaborators, with appropriate recognition given (e.g. co-authorship). Where practical, additional assistance should be solicited ahead of the user's experiment.



Figure 3.1. (A) Overview of the SM Beamline (B) M1 Mirror inside the SM Hutch, (C) EPU inside the Storage ring and (D) SM Beamline Computer and Work Area.



Figure 3.2. STXM Hutch. (A) STXM End Station and Computer, and (B) STXM Control





Figure 3.3 (A) Optical Microscope and computer and (B) User work area.

3.0 Overview of the SM Beamline and STXM Microscope

The SM beamline consists of an elliptically polarizing undulator (EPU) located inside the storage ring, the M1 mirror inside the SM hutch, and a plane grating monochromator (PGM) (Fig.3.1). There are two branches on the SM beamline with dedicated end stations, the scanning transmission X-ray microscope (STXM) (Fig. 3.2) and photoemission electron microscope (PEEM) (Fig. 3.1). The STXM branch has 1 mirror (M3STXM) and the PEEM branch has 2 mirrors (M3PEEM, M4PEEM) (Fig. 3.1). Note that only 1 branch is operational at any one time since the M3PEEM mirror must move in front of the M3STXM mirror in order for the PEEM end station to receive beam.

There are 2 separate SM Beamline Control Programs used to control the beamline for the STXM endstation, namely the SM Beamline Control Graphical User Interface (not shown) and the SM Beamline STXM Branch Graphical User Interface (STXM_Branch_GUI) (see Fig. 4.1). The User is not permitted to use the SM Beamline Control Graphical User Interface unless otherwise directed by the SM Beamline staff.

The STXM end station consists of the STXM Tank (Fig. 3.2A) and the STXM work area (Fig. 3.2B), including the STXM CONTROL Computer and display #3 for the SM Beamline Control computer. The SM Beamline work area is located outside the STXM Hutch, including the SM BEAMLINE CONTROL Computer (the 2 computer displays on the left) (Fig. 3.1D).

Also at the SM Beamline is an Optical Microscope (VLM) for the previewing of samples (Fig. 3.3).

Table 3.1 summarizes the important parameters of the SM beamline. The energy range of the SM beamline is from 130 to 2500 eV (soft X-ray). A table of the elements accessible in the energy range of the SM beamline, as well as the specific edge, energy range, grating, harmonic and polarization for a particular element is shown in Appendix A (The table may not be complete, so if you do not see the element or edge of interest, contact the Beamline Staff for confirmation that the element cannot be run at the SM beamline).

	seumine important parameters.
Source	Apple II type Elliptically Polarizing Undulator, period=75mm
End Stations	STXM, PEEM
Energy Range	130 to 2500 eV
Flux	STXM: 10 ⁸ ph/s for a resolving power of R~3000
	PEEM:10 ¹² ph/s for a resolving power of R~3000
Spatial Resolution	STXM: 30 nm; PEEM: 50 nm
Resolving Power ($E/\Delta E$)	Nominal 3000, can reach $>10^4$

 Table 3.1. Summary of SM beamline important parameters.



Figure 4.1. STXM Branch Control Graphical User Interface (STXM_Branch_GUI) accessed from the SM Beamline Control computer. Whenever possible the STXM_Branch_GUI will be put on the back of the page in Section A.



Figure 4.2. STXM, PEEM, SM, OTHER Button on the task bar located on SM Beamline Control Display 2.

4.0 PART A: Enabling the SM Beamline

The SM Beamline Computer controls the beamline, including the values, shutters, EPU, grating and mirrors. The values and shutters need to be opened before setting-up and optimization of the STXM Microscope.

The following sections show in detail how to enable the beamline and is summarized in a Quick Guide (Section 4.6). The sections are presented in the recommended order for setting up the beamline.

4.1 Starting the STXM Branch Graphical User Interface

Usually, the STXM_BRANCH_GUI (Fig. 4.1) is already running and is visible on the SM Beamline Control Display #3 located in the STXM hutch, right of the STXM Control Monitor. If the STXM_Branch_GUI is already running proceed to Section 4.2, otherwise,

the STXM_Branch_GUI can be accessed from the SM beamline Control Computer, located outside of the STXM hutch (Fig. 3.1D):

- 1. Make sure you are on the STXM display window and not the SM, PEEM or OTHER display window (Fig. 4.2)- button located on the task bar, usually displayed on SM Beamline Control Computer Display #2. This could be the problem if you are not seeing the SM_STXM_GUI.
- 2. To access the SM_STXM_GUI click on the runSTXM shortcut on the desktop display (Display 2). Drag the SM_STXM_GUI onto the SM Beamline Control Computer Display #2, which will also be displayed on the SM Beamline Control Computer Display #3 in the STXM hutch, next to the STXM Control computer monitor.

4.2 Branch Selection

Check to determine whether the STXM or PEEM branch is selected which is indicated in the Branch&Grating Selection section of the SM_STXM_GUI (Fig. 4.1). Change from the PEEM branch to the STXM branch by clicking on the PEEM icon and dragging down to the STXM icon. The changeover is complete when the yellow "not in position" disappears and is replaced by "Branch 2: STXM".

4.3 Safety and Photon Shutters/Access Control Information System Panel

To protect the SM beamline it is divided into sections, separated by valves that can isolate each section automatically in the case of a sudden pressure increase or are used to isolate sections during scheduled maintenance/beamline development. Shutters are located in various spots on the beamline to block X-rays from going further downstream in the beamline.



Figure 4.3. Safety (SSH) and Photon (PSH) Shutters/Access Control Information System Panel. The panel is located at the start of the SM beamline.



Figure 4.1 STXM_BRANCH_GUI

The sample can be loaded prior to opening the shutters and enabling the beamline and/or STXM microscope. The SAMPLE LOADING/UNLOADING procedure is shown in Appendix B.

The Safety (SSH) and Photon (PSH) shutters are closed by the Control Room during an injection, when there is a beam dump or for planned maintenance. These shutters must be opened before beam can come down the beamline. The SSH indicator light/button and the PSH indicator light are located on the Access Control Information System (ACIS) Panel (Fig. 4.3). Permission to open the SSH shutter is given by the Control Room and is indicated when all the lights before the SSH shutter are green. Open the SSH by pushing the LARGE GREEN BUTTON, at which time the LARGE RED BUTTON will go off and the GREEN BUTTON on.

The PSH shutter can only be opened after the SSH shutter is opened, and is opened from the SM-STXM-Branch_GUI (Fig. 4.1). Note that the SSH and PSH shutters also control the beam access to the REIXS beamline and may be opened by their Users/Staff.

If you cannot open the SSH or PSH shutters after the injection is complete, contact the floor coordinator for assistance.

Note also, that the SSH and PSH shutters cannot be closed at the Beamline once they are open.

Figure	4.1	STXM	Branch	GUI
.				

SM Beamline Control System							
10ID-1 SM Beamline Control ST	RingStatus: 0.001 mA Spring Shutdown	2.9 GeV					
PSH1 F.Val. PSH SSH F.Val. VA 1 2 3 4 5 6 7	M1 SM-PSH 4J1 F 8 9 10 1	PGM M3PEEM M3STXM STXM-ES	1 6				
Branch&Grating Selection Note: To change branch or grating, the SM photon shutter (SM-PSH) has to be closed. Peter step#2 SelectBranch PEEM BRANCH 2: STXM SelectGrating LEG MEG (To charge grating, first clock HOME, when it is dore, then select grating) OpenBeamling SSH-CLOSED OPEN/CLOSE	S EPU/Monochromator EPU Gap 200.0000 40.2995 Q1 0.0000 0.0010 Q2 0.0000 -0.0025 Q3 0.0000 -0.0025 Q4 0.0000 -0.0010 Taper 0.0000 0.0115 MONOCHROMATOR Mirror 36415.65 36419.39 Grating 49705.86 49705.99 • 6 REIXS EPU	8 Valves&Vacuum Sec.1 07 1.166e-10 Sec.2 0P 9.829e-11 Sec.4 07 01 1.367e-09 Sec.6 0P 01 1.367e-09 Sec.6 0P 01 6.874e-10 Sec.7 0P 01 4.807e-10 Sec.8 0P 01 3.885e-10 Sec.9 0P 01 6.148e-10 Sec.10 0P 01 6.148e-10 Sec.11 0P 01 6.761e-10 Sec.12 0P 01 5.026e-10 Sec.13 0P 01 1.089e-10 StXM Sec.14 0P 01 1.095e-09					
If you want the EPU or grating tracking OFF during a scan, Please contact the beamline staff. Hammonic Mode Angle 1 LintHor 0.000 0.000 Cff 2.149 Epu-offset 0.000 Energy (eV) 320.052	Gap 41.2145 Q1 0.0035 Q2 -0.0025 Q3 0.0020 Q4 -0.0020 Taper 0.0345 7 UsefulTools Scan Picoammeter	Sec.15 OP CI 3.970e-10 •• Sec.16 OP CI 5.210e-10 •• Pressure should be < 5E-9 before opening any valve					

4.4 Beamline Safety and SM-Photon Shutters and the Variable Aperture

The SM-PSH shutter is to isolate and protect the PGM and the beamline safety shutter (BSH) is to isolate the STXM Tank and protect the User from the X-rays.

The SM-PSH shutter is to be closed whenever there is a grating or branch change.

The BSH shutter is to be closed whenever there is a sample change or the STXM Tank is being accessed.

1. To Open/Close the SM-PSH or BSH shutters click on the OPEN or CLOSED buttons.

Also, there is a variable aperture (VA) before the M1 mirror that protects the mirror from the beam when the beamline is not in use. The variable aperture closes when the CloseBM button was clicked or can also be closed by the REIXS Users/Staff. However, it is not apparent from STXM_Branch_GUI valve display that the variable aperture is closed- is apparent on the SM-BEAMLINE_GUI valve display (not shown) as the VA valve is RED when closed.

2. To open the variable aperture (VA) click on the ``OpenBM`` button.

There should now be beam to the STXM End Station. Follow the procedures in Section 5 to verify that there is beam to the STXM Microscope.

4.5 Closing the SM Beamline (at the end of the shift)

At the end of the shift we want to prevent beam from hitting the mirrors, particularly the M1 and silicon nitride window on the zone plate snout, as this will extend their life and cut down on carbon contamination. This is accomplished by closing the shutters and variable aperture.

- 1. Close the SM-PSH shutter
- 2. Close the BSH shutter
- 3. Close the variable aperture by clicking on the "CloseBM" button.
- 4. Turn OFF the detector.

4.6 Enabling the SM Beamline (Quick Guide)

1	
	Enabling the SM Beamline (Quick Guide)
	(STXM_BRANCH_GUI, unless otherwise stated)
	4.0 Enabling the SM Beamline
	(Sample can be loaded at anytime)
	4.1 Start STXM Branch GUI- usually already displayed
	4.2 Branch Selection (Close SM-PSH shutter)- Select STXM branch if on PEEM Branch
	4.3 Open safety shutter (only after injection) – Green button on ACIS panel and Open PSH shutter
	4.4 Open SM-PSH and BSH shutters and Click on OPENBM (opens variable shutter)
	4.2 Grating or Branch Change
	1. Close SM-PSH shutter
	Branch change (only necessary if on PEEM branch)
	a. Click PEEM and drag down to STXM
	Grating change
	a. Click HOME– when complete indicates DONE
	b. Click the appropriate grating (e.g., LEG, MEG, HEG) – when
	complete a number is indicated (250=LEG, 500=MEG, 1250=HEG)
	2. Open SM-PSH shutter

Sample Change (Appendix B)

- 1. Close SM-PSH and BSH shutters
- 2. Turn OFF the Detector
- 3. Move "Coarse Z" to 5000 um
- 4. Open valve to let in air into the STXM End Station

4.5 Closing the Beamline – at the end of the shift

- 1. Close SM-PSH and BSH shutters
- 2. CloseBM (closes variable aperture)
- 3. Turn OFF the Detector



Figure 5.1 STXM Control window when the computers restarts.

Table 5.1. Preferred grating and harmonic settings as a function of photon energy.

Energy Range (eV)	Grating	Harmonic	Polarization
130 to 390	LEG	1^{st}	LinHor, CirRight, CirLeft
330 to 1000	MEG	1 st	LinHor, CirRight, CirLeft
1000 to 1800	MEG	3 rd	LinHor
1800 to 2500	MEG	5 th	LinHor
>2200	HEG	5 th or 7 th (not well characterized)	LinHor



Figure 5.2 Photon flux as a function of energy at the STXM exit slits for (A) First order light using circular polarization, and (B) First, third and fifth order light using linear horizontal polarization.

5.0 PART B: STXM Microscope Set-up

The following assumes that the SM beamline has been enabled.

The STXM computer controls the microscope set-up and optimization and data acquisition (See STXM User Manual for details on data acquisition), and is located inside the STXM hutch (Fig. 3.2A).

Sometimes it is necessary to restart the computer, which is done either through the Windows software or by pushing the button located on the front of the computer (Fig. 3.2A)

If the STXM computer is restarted, the User must log back on and the screen opens to Figure Fig. 5.1. The login selection and password (PW) are found on the bottom left of the STXM Control computer screen monitor, and the PW is capital sensitive.

Note, usually the sample is loaded previous to setting up and optimizing the STXM Microscope. To see the loading procedure, as well as the sample plate design and numbering system see Appendix B.

Also, Users need to understand the relationship between the zone plate, order sorting aperture (OSA), Sample and Detector, which is discussed in Appendix D. We encourage new Users to look at Appendix D before proceeding.

5.1 Zone Plate Selection

There are a number of Zone Plates that can be used, defined by their spatial resolution, mainly 25, 35 and 40 nm. For energies before ~280 eV the 40 nm zone plate must be used as the focal length of the 25 and 35 nm prevent their use due to physical constraints. Only Beamline Staff can change zone plates. It takes about 30 mins to change a zone plate and set-up the microscope again. Experiments should be planned so that the changing of zone plates is kept to a minimum.

5.2 Grating Selection

The SM beamline has 3 gratings, the low energy grating (LEG, 250 l/mm), the medium energy grating (MEG, 500 l/mm) and the high energy grating (HEG, 1250 l/mm). The flux of the gratings is energy dependent, thus, choosing the appropriate grating depends on the element of interest. Table 5.1 shows the photon energy range of each grating and the photon flux as a function of energy is shown in Figure 5.2.

Figure 4.1 STXM_Branch_GUI

SM Beamline Control System					
10ID-1 SM Beamline Control ST	XM branch	RingStatus: 0.001 mA Spring Shutdown	2.9 GeV		
PSH1 Val. PSH SSH F.Val. VA 1 2 3 4 5 6 7	M1 SM-PSH 4J1	PGM M3PEEM M3STXM STXM-ES	16		
Branch&GratingSelection Note: To change branch or grating, the SM photon shutter (SM-PSH) has to be closed. Refer step#2 SelectBranch PEEM O BRANCH 2: STXM ScleatGrating UD	5 EPU/Monochromator EPU Gap Gap 200.0000 40.2995 0 Q1 0.0000	8 Valves&Vacuum Sec.1 00 1.166e-10 • Sec.2 00 9.829e-11 • Sec.4 00 CL 1.367e-09			
Cocharge grating test dick Holding HEG 250 (To charge grating thist dick Holding) when it is done, then select grating) 2 OpenBeamline SSLI CLOSED PSH-CLOSED	Q2 j0.0000 -0.0025 Q3 j0.0000 0.0025	Sec.6 OP CL 6.674e-10 ••• Sec.7 OP CL 4.807e-10 Sec.8 OP CL 3.885e-10 •••			
OPEN/CLOSE OPEN CLOSE OPEN CLOSE	Q4 0.0000 -0.0010 O Taper 0.0000 0.0115	Sec.9 OP CL 5.148e-10 000			
3 CloseBeamline	Minror 36415.65 36419.39 Grating 49705.86 49705.99	Sec.12 OP CL 5.026e-10 000 Sec.13 OP CL 1.089e-10 000 STXM			
If you want the EPU or grating tracking OFF during a scan, Please contact the beamline staff.	Gap 41.2145 Q1 0.0035 Q2 -0.0025	Sec.14 OP CL 1.095e-09 000 Sec.15 OP CL 3.970e-10 00 Sec.16 OP CL 5.210e-10 000			
Hacmonic Mode Angle LimHor [0.000] 0.000 Cff 2.149 Epu-offset 0.000	Q3 0.0020 Q4 -0.0020 Taper 0.0345	Pressure should be < 5E-9 before opening any valve Updated: 26Jan2011 Contacts:			
Energy (eV) 320.052	7 UsefulTools Scan Picoammeter StripTool	Jian Wang: 3546 Yingshen Lu: 3743 Jay Dynes: 3840			

Note, changing the grating will result in the energy being changed.

- 1. To change gratings first click on the "HOME" button on the STXM_Branch_GUI (Fig. 4.1), which moves the PGM tank to a set position. The display will flash "not in position" while the tank is moving, and "DONE" when the tank has stopped moving. Note it may be necessary to click the ``HOME`` again as sometimes the move is not completed as evident by `DONE` not appearing and the `not in position`` stops flashing.
- Then click on the desired grating (i.e., LEG, MEG or HEG) (Fig. 4.1), and while the grating is moving the display will show "not in position" and will display the numbers 250 (LEG), 500 (MEG) or 1250 (HEG) when complete.

Following this procedure ensures that the energy calibration is more reproducible. If the grating is selected without pushing the "HOME" button the energy scale may change by up to 0.3 eV from its previous setting when changing gratings.



Figure 5.3. STXM Control Graphical User Interface (STXM_CONTROL_GUI). Many of the menus have pull downs. Whenever possible the STXM_CONTROL_GUI will be put on the back of the page in Section B.

TOG	us & zone Plates Parameters						-	
	Focal Length Function A1 A0 -4.86576 400	-Zone Plat	Diameter	Central stop	Outer zone	A1		Current Zone Plate
	f(E) = -4.866E + 400	• # 1 • # 2 • # 3	240	90	35	-6.7916	••	
ŝ	A0 is OSA to SAMPLE distance (um)	C # 4	150	55	60	-7.2692		
	Expected positions for current energy Focal Length (µm) 1557.29	OSA • #1 • #2	Diameter 50				J	
	Sample Z (add sample thickness) (μm) 400 Zone Plate Z (μm) -1157.29	C #3 Max Ao for	60 current ene	ergy estimation	ation: Ao <	308		
	Calculated Zone Plate	Positi	on		1	OK	1	

Figure 5.4. Zone Plate focus parameters and Ao (OSA to Sample distance (um)).

5.3 Energy Selection

<u>CAUTION:</u> There is a risk of crashing the zone plate into the Order Sorting Aperture (OSA) when going from HIGHER ENERGY to LOWER ENERGY. The next couple of steps are to protect the zone plate from crashing into the OSA.

We will now use the STXM Control Graphical User Interface (STXM_CONTROL_GUI) (Figure 5.3) to change the energy.

- 3. Change the zone plate `focus` from AUTO to STATIC (``focus`` window, bottom middle) (Fig. 5.3). This is done by dragging from the AUTO to STATIC. In the static mode, the zone plate will not move when a new energy value is entered.
- 4. Enter 395 eV in the Destination box (or the Energy for your experiment) and either press the ``Go`` button or Enter. This will calculate the position that the zone plate should move to using the current Ao value. Note that the Ao value is the expected sample position relative to the OSA. Thus, the focal length (f) is calculated based on the distance of the zone plate from the OSA (Z) and the distance from the OSA to the Sample (Ao) (f = Z + Ao). Ao increases with increasing energy. The reason that we change the Ao instead of just moving the zone plate further back from the OSA is due to alignment issues and physical constraints. So if we move from higher energy to lower energy, the calculation of the new focal length is based on the Ao from the higher energy, subsequently the zone plate may crash into the OSA when moving from higher energy to lower energy.
- 5. Change Ao to the recommended value. To see the recommended Ao value, click on the ``S`` box next to the ``Focus`` box (AUTO or STATIC displayed) (Fig. 5.3), which opens into a new window (Figure 5.4). The maximum Ao for the current energy estimation is shown on the bottom right. Select an Ao value less than the maximum, usually within about 100 um of the Ao maximum. The exact Ao value is not that critical as the zone plate position will adjust accordingly. Input the new Ao (top left) and click OK. Also, notice the zone plate that is currently being used, indicated by the round circle being full.
- 6. Change the zone plate focus back from STATIC to AUTO. If this is not done, the zone plate position will not change with a change in energy.
- 7. Click on the energy `Go`` button. Observe that the Zone Plate position changes. At this time it would be advantageous to look through the STXM Tank viewport and estimate the distance between the zone plate and the OAS. The expected value is that shown by the ZonePlate box. If you think the zone plate and OSA are touching you may need to open the STXM Tank. Note, if the zone plate and OSA are touching, and you continue to acquire data, irreversible damage may result to the zone plate and/or the OSA. Call the Beamline Staff to assist you before damaging that \$10,000 zone plate.

_ 🗆 🗙 orage Ring Current * A110115007.hdr Scan mA On Off 🖬 🎒 🗱 🕄 Set as lo counter0 mage Scan 200x200 pts 215.004 eV 2.095 ms Dwell (f) counter0 -111000 30 \$ 1055 110000 108000 25 106000 20 104000 15 102000 100000 98000 96000 38 95000 293.5 294 0 • Data 294.029 95506.6 Sample X 0 0 Abs Energy counter0 Sample Y No Data Sample To Cursor Loading scan: A110115022.hdr Reading file A110115022_0.xsp E:\Data\STXM-data\2011\2011_01\110115\A110115022.hdr loaded success Loading scan: A110115007.hdr Reading file A110115007_0.xsp E:\Data\STXM-data\2011\2011_01\110115\A110115007.hdr loaded success 1055 ٠ PMT /10 Sample To Cursor - 800 Coarse - 600 Auto - 400 7 138 30 Updates On Microscope Status SampleX • SampleY • SampleFineZ • nline Control 0.000 5000.395 -25.00 0.00 5000.7 5100.1 -968.7 6 Destination Current Scan Status Scan Control Energy @ 1541.720 @ 000000000 1473 Go S of No Current Scan Image -(eV) START Sample Scan Estimated Time of Region Hori_Defl_(urad) • -67.6 • -67.6 Go S CoarseX CoarseY CoarseZ of Elapsed Time Line r Exit Slits -Dispersive (µ Next Scan # Point of 0 38 38 0 GO S ZonePlateZ 111 0 111 Microscope Control Motor OSAX ... NonDispersive(µ 0 Go S Current Destinati Setup Motor n / Jog OSAY 00000000 All OFF All ON 0.5 -0.6 11.7 100.4 0.8 1 0 60 - 0 -0 Go D Stop S CoarseX 0.00 • DetectorX DetectorY 0 69 Polarization • • 5100.1 @ 5100.1 Go D Stop S CoarseZ DetectorZ Cff det Off EPU Folowing Con EPU Harmonic 1st Focus -
 Sample
 OSA IN
 Zone Plate IN
 Sample IN

 Mover
 OSA OUT
 Zone Plate OUT
 Sample OUT
 Scan Crd Auto • S Reset Sto **BL** Feedback Grating Interfer AUX3 AUX4 0.0 0 Exit Slit Cur it Curr Ao 400.0 Acquire lock 0.000 Fbk Offse 0.000 0.000 set 0.000 0.000

Figure 5.3 STXM_CONTROL_GUI

Codes for Polarization Setup for STXM

Table 5.2. EPU – Vertical, Circular, and Horizontal Polarizations

Actual Polarization	IncVert (-90)	CirLeft	LinHor	CirRight	IncVert (+90)
STXM value	-2	-1	0	1	+2

Table 5.3. EPU-Inclined Polarizations¹

Actual	-80	-70	-60	-50	-40	-30	-20	-10	10	20	30	40	50	60	70	80
Angle																
STXM	20	30	40	50	60	70	80	90	110	120	130	140	150	160	170	180
value																

1. Add 100 to the actual angle required, -80 (Actual Angle) + 100 = 20 (STXM value entered).

5.4 Harmonic Selection

The harmonic determines the EPU gap settings, which influences the flux. The recommended harmonic for each energy range is displayed in Table 4.1 or for each element in Appendix A.

1. To change the harmonic drag down to the desired harmonic using the EPU Harmonic box on the STXM_Control_GUI (Fig. 5.3) or click on the HARMONIC button on the STXM_Branch_GUI (Fig. 4.1). It is preferable to use the STXM_Control_GUI as the parameters file updated, whereas making the change in the STXM_Branch_GUI does not update the parameter file.

Note, initial optimization of the microscope is usually done at 395 eV (just below the N K-edge), so initially use the values for N K-edge (i.e., harmonic 1).

5.5 Polarization Selection

The SM beamline is designed to produce plane horizontal (LinHor) or vertical (IntVert- (90°) , IntVert+ $(+90^{\circ})$, circular left (CirLeft) or right (CirRight) and linear polarized light (LinInc) with the polarization planes inclined at any angle (see Appendix C for more details). For nonmagnetic samples, where no change in polarization is required, 3 settings are available (LinHor, CirLeft or CirRight). The polarization influences the flux at a selected energy, thus, the polarization selection depends on the element of interest. The recommended polarization for each energy range is shown in Table 5.1 or for each element in Appendix A.

Note, initial optimization of the microscope is usually done at 395 eV (just below the N K-edge), so initially use the values for N K-edge (i.e., polarization either CirLeft or CirRight)

1a.To change polarization enter the appropriate value in the EPU Polarization box (CirLeft=-1, CirRight=1, LinHor=0, IntVert -90°=-2 and IntVert +90°=+2) (Table 5.2) on the STXM_Control_GUI (see Fig. 5.3) or click on the MODE button in the STXM_BRANCH_GUI, and drag down to the appropriate polarization. The advantage of making the change in the STXM_Control_GUI is that the polarization used is recorded with the data, otherwise the previous polarization is recorded with the data.

For magnetic samples, where a change in polarization is required (e.g., X-ray magnetic circular dichroism (XMCD) studies), the angle can be varied from 20 to 180 degrees.

1b. To change polarization enter the appropriate value in the EPU Polarization box (Table 5.3) on the STXM-Control_GUI or click on the MODE button, and drag down to LinInc and enter the appropriate value for the angle. To obtain the actual angle required a value needs to be entered into the ANGLE box. The value to enter is determined by adding 100 to the angle required. For example, for an actual angle of -80° , enter the value 20 (-80 + 100). See Table 4.3 for the values for angles for every 10° , starting at -80° .

5.6 Constant of Fixed Focus (Cff)

The constant of fixed focus (Cff) value controls the angles for the M1 mirror and the grating of the PGM, which determines the amount of the reflected and diffracted light hitting the M3STXM mirror. Generally the higher the reflected light, the better the spectral resolution, while the higher the diffracted lighted the higher the flux. Usually the Cff value is set at 2.15, a value that gives a resolving power over 3000 and good flux. By increasing the Cff value, the resolving power will increasing at the expense of the flux (Appendix H). The Cff value is displayed on the STXM_BRANCH_GUI but can only be changed from the SM_CONTROL_GUI._To change the Cff value consult with Beamline Staff.



Figure 5.3. STXM_CONTROL_GUI



Figure 6.1 Photomultiplier Tube (PMT) High Voltage Controller on top of the STXM Tank. Remember to turn off the High Voltage when opening the STXM Tank to light.

6.0 PART C: STXM Microscope Optimization

To optimize the STXM Microscope requires that the flux be maximized for a particular set-up.

This requires that some beam reach the detector.

6.1 Verify that Beam is reaching the Detector

It is assumed that the beamline has been enabled, that is, all valves and shutters are open including the variable aperture.

The STXM microscope is routinely optimized at 395 eV as there is lots of flux and the amount of flux expected is well known. Thus, to verify that the detector is responding to the beam, we will use the set-up conditions for 395 eV (See Section 5 for details).

Grating: MEG Harmonic: 1 Polarization: CirLeft or CirRight Energy: 395 eV Exit Slits: 20/20 (Dispersive/Nondispersive) (see below for details)

Also, it is assumed that a sample plate with samples has been loaded (See Appendix B for the procedure for Loading/Unloading a sample) and that the STXM tank has been pumped and filled with 1/6 atmosphere of He.

- 1. Turn on the PMT detector by turning on the high voltage controller, which is usually sitting on top of the STXM chamber (Fig. 6.1). The voltage should be around 1050 V.
- 2. Go to an open hole on your sample plate (assumes a sample is loaded). See Appendix B for the sample plate coordination system.
- 3. Move the sample holder to 3000 um using ``Coarse Z``. This distance is just to ensure that the sample holder is not blocking the beam. For looking at actual samples, the expected position between the sample holder and OSA is about +300 um larger than the Ao value (e.g., Ao = 300, expect Coarse Z about +600). CAUTION is advised though, as sometimes the calibration of the Coarse Z scale goes off. Thus, it is usually better to gradually move (e.g., starting position = Ao + 1000 um) the Coarse Z closer (i.e., sample holder) and use the focus as the guide.
- 4. On the spectral display window (top right corner) on the STXM_CONTROL_GUI the counts from the detector are displayed. If you see a spectrum on a yellow background it means that the detector counts are not being displayed. Click on the CHART button (very top right) to



Figure 5.3 STXM-Control_GUI



Figure 6.2. (A) Detector scan parameters and (B) detector scan image

- get the count screen. If the detector is connected to the STXM computer there should be counts around 400 in the display, even if the SM-PSH and BSH shutters are closed. If the counts are reading 0 (i.e., flat line as shown) it means that the detector is not connected to the STXM computer. Consult with the Beamline Staff for possible solutions.
- 5. Open the shutter on the STXM_CONTROL_GUI by changing the ``AUTO`` to ``OPEN. Note, the shutter is usually left in the ``AUTO`` `position as it is closed in the AUTO setting, opening automatically when acquiring data. There should be an increase in the number of counts from when the shutter was closed. Usually we expect 200,000+ counts at 395 eV, depends on exit slit size and whether the microscope has been optimized. Nevertheless, as long as the counts were observed to increase when the shutter was opened (e.g., 2000) it should be possible to optimize the microscope using the procedures outlined in the next section.
- 6. If there was no increase in the counts, it indicates that the beam is possibly being blocked. To trouble shoot consider the following:
 - a. PMT detector is turned off;
 - b. You are not in an open hole on the sample plate;
 - c. Valves are closed;
 - d. OSA is not aligned (see below);
 - e. Detector is not aligned (see below).

6.2 Detector Scan

For the best sensitivity, the detector needs to be very well aligned with the beam.

The Detector Scan menu can be accessed from the Scan Control by dragging down to the Detector Scan (Fig. 5.3). The usual parameters for a detector scan are shown in Figure 6.2A. Click "Begin Scan" and the expected result of a well centered detector is shown in Figure 6.2B. To align the detector, move the cursor to the centre of the Detector image by clicking on the detector image, which will activate the "Detector to Cursor + Set to 0,0" button on the left hand side under the image display window. Click on the "Detector to Cursor + Set to 0,0" button.

If the detector is not found there are likely three reasons. First, make sure you are in a hole. Secondly, the OSA may not be aligned and is blocking the beam. Move the OSA out of way by clicking on the "OSA Out" button. Do another detector scan. Thirdly, increase the scan size from the usual 1000 x 1000 um to 2500×2500 um. Do another detector scan. If the detector is still not found contact the Beamline staff for assistance.

Note that the detector scan usual only needs to be done once, at the beginning of the shift, to confirm that the detector is well aligned.

05A Scan	- Estimated Time:	47s						×
A	Centre Pos (µm)	Range (µm)	# Points	Step (µm)	- inergy (eV) 320.051 ZP in Focus	Begin	Save Scan	Recorded
OSA X	0	70	30	2.333	Dwell Time (ms) 3	Joan		chameisii
OSA Y	0	70	30	2.333	ZP to in-focus	Cancel	Load Scan Definition	
	1							
B				C	D			

Figure 6.3. (A) OSA scan parameters, (B) OSA scan at 395 eV (C) OSA scan at 395 when ZP in focus (is checked) and (D) OSA scan that shows the OSA and Zone Plate (300 x 300 um).

OSA Centre F	'os (μm) Le	ength (μm) Theta (*)	# Points	Step (µm)	Energy (eV)	320.052	Begin Scan
X -30.0 Y 0.0	000	50.0	30	1.667	Dwell (ms)	3	Cancel
Zone Plate	Pos (um)	Range (um)	# Points	Step (um)			Save Scan Definition
Z -15	57.3	200.0	30	6.66667		Recorded Channels	Load Scan Definition
1						<u></u>	r
	В						C
							C

Figure 6.4. (A) Normal OSA focus scan parameters, (B) Green line on OSA image (C) OSA focus scan

6.3 OSA Scan

The purpose of the OSA is to filter out zero order light (Appendix D). The OSA needs to be well aligned with the beam to be most efficient.

The OSA Scan menu can be accessed from the Scan Control by dragging down to the OSA Scan (Fig. 5.3). The usual parameters for an OSA scan are shown in Figure 6.3A. Make sure that the ``ZP in focus`` box is not checked at this time. Click "Begin Scan" and the expected result of a well centered OSA at 395 eV is shown in Figure 6.3B. To align the OSA, move the cursor to the centre of the OSA Scan image by clicking on the OSA image, which will activate the "OSA to Cursor + Set to 0,0" button on the left hand side under the image display window. Click on the "OSA to Cursor + Set to 0,0" button.

If the OSA is not found with the usual 70 x 70 um scan size, the likely reason is that it was bumped during sample loading, increase the size to 300 x 300 um and do another OSA Scan. Centre the OSA as outlined above. The accuracy of centering the OSA is better at 70 x 70 um than 300 x 00 um, so repeat at 70 x 70 um. If the OSA is still not found contact the Beamline staff for assistance.

To check the size of the OSA hole, check the `ZP in focus`` box and the expected result is shown in Figure 6.3C. The size of the OSA can be determined by clicking on the left and right OSA Focus Scan edges of the OSA scan. Generally, the 50 um OSA hole is used. DO NOT Align the image- just for reference.

To see the OSA and Zone plate the range is increased to 300 x 300 um (Figure 6.3A). DO NOT Align on 300 x 300 image- just for reference/ uniformity of the zone plate.

Note, the OSA scan should be done whenever loading a sample to verify that it was not bumped changing the sample.

6.4 OSA Focusing Scan

The OSA Focus Scan is done to calibrate the Zone Plate scale using the OSA as the reference point.

The OSA Focus Scan is done right after the OSA Scan.

The OSA Focus Scan menu can be accessed from the Scan Control by dragging down to the OSA Focus Scan (Fig. 5.2). The usual parameters for an OSA focus scan are shown in Figure 6.4A. Note a Green Line came onto the OSA image when the OSA focus scan was opened (Figure 6.4B). Usually it is not necessary to move the Green line. Click `Begin Scan`` and expected image of a calibrated zone plate is shown in Figure 6.4C. Move the cursor to the focused position by clicking on the OSA focus scan image, which will activate the "Focus to

10ID-1 CLS SM Beamline Manual

Version: 13 June 2011

Figure 5.3 STXM_CONTROL_GUI



Cursor + Set ZP calibration`` button, again on the left hand side under the display window. Click on the "Focus to Cursor + Set ZP calibration``

If it was not possible to see the focal point, it may be because the zone plate was way out of focus. Consult the Beamline Staff for assistance.

The OSA Focus Scan usually needs to be conducted once at the beginning of a shift.

6.5 Exit Slits

The exit slit controls the size of the beam, in both the vertical and horizontal directions, that hits the zone plate. The size of the beam determines both the flux and the spectral resolution. Generally, the larger the beam the more flux but the poorer the spectral resolution. The vertical slit (Dispersive) influences the spectral resolution to a greater extent than the horizontal slit (Nondispersive). There is always a trade-off between the flux and spectral resolution. However, the exit slit size is usually selected for the maximum flux (< 20 MHz, see optimizing section).

The exit slit size is often displayed as the Dispersive/Nondispersive slits (e.g., 20/20).

The size of the exit slits are adjusted to prevent saturation of the PMT detector as a linear response is necessary for quantification and good image quality. The PMT saturates when the Io signal intensity is above 20 MHz/s. Click on the CHART button (very top right) to get the count screen. The chart units are Hz/s. The PMT is usually set to **PMT/10**. So a chart reading of 2 x 10^6 Hz/s is equal to 20 MHz (multiply by 10 for the PMT and 1 MHz = 10^6 Hz). On your image, for 2000 Hz/ms maximum = 20 MHz/s for a dwell time of 1 ms (2000 Hz/ms x 10 (PMT/10) x 1000 ms/s). So when counts > 2000 on an image may indicate that your detector is saturated (if your frame is in the image counts may be > 2000).

Since the signal intensity is a function of the energy, the selection of the energy at which the exit slit size is adjusted to keep the intensity < 20 MHz/s is usually done at the lowest energy of the edge that you will collect a stack/line scan, as generally the Io signal intensity decreases with increasing energy (Figure 4.4). For example at the Fe L-edge, use 695-700 eV. For the C K-edge we generally set-up at 300 eV because of the big dip in the C K-edge Io signal due to C contamination on the optics and windows.

6.6 M3STXM Mirror Pitch

The position of the beam on the M3STXM mirror changes slightly with a change in energy. That is, the beam moves off of the best spot on the M3STXM mirror, resulting in a reduction in the intensity of the beam. The lose in intensity can be corrected for by adjusting the pitch of the M3STXM Mirror.

Figure 5.3 STXM_CONTROL_GUI



Figure 6.5. (A) Hori_Defl_(urad) (M3STXM pitch) window open by clicking on "S"

The M3STXM pitch is shown on the STXM_CONTROL_GUI as the Hori_Defl_(mrad) box (Fig. XX). Click on the "S" button beside the Hori_Defl box, opening to Figure 6.5. Under the JOG, enter 1 if not already entered (sometimes 0 if STXM_Control program restarted). Clear the display (Clear button top left of display) to make it easier to see small changes. Use the red arrows left and right of the JOG box to maximize the intensity. If the signal is going down, use the other arrow. Be patient between each movement, as there is a backlash sometimes and the signal may actually decrease initially then increase. When optimized, close the Hori_Defl "S" display. You may have to readjust the exit slits so that the Io signal maximum < 20 MHz/s.

6.7 EPU Offset

The position of the beam relative to the EPU Gap changes slightly with a change in energy. That is, the beam may be off centre of the EPU Gap, resulting in a reduction in the intensity of the beam. The lose in intensity can be corrected for by adjusting the EPU Gap Offset.

To center the beam down the EPU is done by changing the EPU Offset (Fig. 5.3). Adjust the EPU offset to maximize the intensity. Clear the display (Clear button top left of display) to make it easier to see small changes. The EPU offset is usually between -0.4 and 0.4. You must click on the "SET" button to implement the change. Confirm that the change was made by comparing that the actual and input values are the same. Make sure the actual value and input value are the same. Adjust by 0.05 units until the maximum signal is realized. You may have to readjust the exit slits so that the Io signal maximum < 20 MHz/s.

6.8 STXM Microscope Set-up and Optimization Summary

STXM Microscope Set-up and Optimization (Quick Guide)
Use STXM_Control_GUI
Consult the manual for complete details
 STXM Set-Up 5.1 Select Zone Plate (Discuss options with Beamline Staff) 5.2 Select Grating 5.3 Select Energy CAUTION : Changing from High to Low Energy. a) Change zone plate focus from "AUTO" to "STATIC" b) Enter Energy c) Change Ao (click on "S" to open new screen) d) Change zone plate focus from "STATIC" to "AUTO" e) Click Energy "Go" button 5.4 Select Harmonic 5.5 Select Polarization 5.6 Select Cff (usually 2.15)
 STXM Optimization 6.1 Verify the Beam is reaching the Detector Change "AUTO" to "OPEN" to open the shutter If no beam check that: a. PMT detector is turned on; b. You are in an open hole on the sample plate; c. All valves are open; d. OSA is aligned (see below); e. Detector is aligned (see below). 6.2 Detector Scan 6.3 OSA Scan 6.4 OSA Focusing Scan 6.5 Adjust Exit Slits (< 20 MHz/s counts, 2 x 10⁶ on display) 6.6 Adjust M3STXM Mirror Pitch (i.e., Hori_Defl_(uram) 6.7 Adjust EPU Offset



Figure 7.1. (A) Gas Cylinders at the SM entrance, (B) Vacuum pump switch and gas inlet values for He, CO_2 , Ne and N_2 , (C) Analogy Gas gauge for measuring He level in the STXM Tank, (D) Digital gauge for measuring gases pressure for calibration, and (E) Vacuum pump value.

Table 7.1 Gas Calibration Lines for STXM

Edge	Gas	Line	Energy, eV	References
S 2p	SF ₆	A_{1g} (1/2)	173.44	Hudson et al, Phys Rev A 47 (1993) 361
		4s Ryd (3/2)	177.42	
		T_{2g} (1/2)	184.57	
S 1s	SF ₆			
C 1s	CO_2	$\pi^{*}(v=0)$	290.74	Ma et al, Phys Rev A 44 (1991) 1848
		3s(v=0)	292.74	
		3p(v=0)	294.96	
N 1s	N ₂	$\pi^{*}(v=0)$	400.87	Chen et al, Phys Rev A 40 (1989) 6737
		3s(v=0)	406.150	
		3p(v=0)	407.115	
O 1s	O ₂	$\pi^{*}(v=3)$	530.75	Prince et al, J. Elec. Spect. 42 (1999) 141
		Rydberg	541.85	Ma et al, Phys Rev A 44 (1991) 1848
	CO_2	$\pi^*(v=0)$	535.4	Hitchcock et al, J. Elec. Spect. 42 (1987) 11
		3s(Ryd)	538.9	
F 1s	SF ₆	$6a_{1g}$	688.0	Francis et al, Phys Rev A 52 (1995) 4665
		$2t_{2g}$	698.9	
Ne 1s	Ne	3p	867.12	Coreno et al, Phys Rev A 59 (1999) 2494
		4p	868.69	

7.0 SM Beamline (STXM) Energy Calibration

The energy scale of the beamline is calibrated generally about 3 or 4 times a year by the Beamline Staff, taking about 4 hours to complete. However, for many reasons, the energy scale may change after the calibration procedure. Thus, Users are recommended to use reference compounds (with known peak positions) and/or compressed gases to accurately adjust the energy scale of their collected data.

Compressed gases available for calibration at the beamline include carbon dioxide, nitrogen, neon and SF₆. The **carbon dioxide**, **nitrogen**, **and SF₆** are outside the STXM hutch near the entrance to the SM beamline, while Ne is in a small cylinder inside the STXM hutch (Fig. 7.1).

Gas Calibration Procedure (Note, only Beamline Staff are permitted to change the gas cylinders)

- 1. Identify the gas that is closest to your energy of interest (Table 7.1).
- 2. Remove the sample or move to a hole in the sample plate (Gas calibration can be done with your sample inside the STXM Tank).
- 3. Pump the STXM Tank to below 200 mTorr, making sure that the pump value is closed when the pump is first started, then open the value after the pump is started. Note the End Station should be opened to atmosphere if it was filled with He, as He is very difficult to pump.
- 4. Fill the STXM Tank with the recommended pressure of gas (Table 7.1) using the digital gauge. If too much gas is allowed into the Tank, pump the excess gas out until the desired pressure is reached.
- 5. Select the Exit Slits size (Table 7.1).
- 6. The calibration parameters are found in the Folder "Gas-Calibration". Use the STXM_CONTROL_GUI to navigate to them.
- 7. The Gas Calibration Lines for STXM are shown in Table 7.2 and examples of the spectra in Figure 7.2 (next page).

Edge	Gas	Pressure (Torr)	Exit Slit Size (µm)					
Carbon	CO ₂	3	30x30					
Nitrogen	N_2	2	30x30					
Oxygen	O_2 (Air) or CO_2	7	30x30					
Neon	Ne	20^{1}	10x10					
Fluorine	SF ₆	3	15x15					
Sulfur	SF ₆	XX	XX					

Table 7.2. Parameters for the Compressed Gases

1 – The gauge has trouble reading the Neon, so fill to 10 Torr, and it will continue to increase.



Figure 7.2. Gas spectra used for calibrating at the C, N, O, F and Ne 1s edges.

Appendix A: Selecting Parameters for Specific Elements

Element	Edge	Energy	Grating	Harmonic	Polarization	Comments
	C	Range (eV)	C C			
В	Κ	185-220	LEG	1^{st}	LH, CR, CL	
С	Κ	280-320	LEG	1^{st}	LH, CR, CL	CO ₂ for calibration
Ν	Κ	390-460	LEG	1^{st}	LH, CR, CL	N ₂ for calibration
0	Κ	520-560	MEG	1^{st}	LH, CR, CL	Air/CO ₂ for calibration
F	Κ	675-715	MEG	1^{st}	LH, CR, CL	SF ₆ for calibration
Ne	Κ	865-870	MEG	1^{st}	LH, CR, CL	Calibration gas
Na	Κ	1060-1090	MEG	1^{st}	LH	
Mg	Κ	1295-1320	MEG	1^{st}	LH	
Al	Κ	1550-1610	MEG	$3^{\rm rd}$	LH	
Si	Κ	1835-1890	MEG	5 th	LH	
Р	Κ	2135-2165	MEG	5 th	LH	
Р	L	125-155	LEG	1^{st}	LH, CR, CL	
S	Κ	2415-2500	MEG			SF ₆ for calibration
S	L	155-190	LEG	1^{st}	LH, CR, CL	SF ₆ for calibration
Cl	L	190-215	LEG	1^{st}	LH, CR, CL	
Ar						
Κ	K	290-320	LEG	1^{st}	LH, CR, CL	
Ca	L	340-360	LEG	1^{st}	LH, CR, CL	
Sc						
Ti	L	450-490	MEG	1^{st}	LH, CR, CL	
V	L	505-560	MEG	1^{st}	LH, CR, CL	
Cr	L	570-600	MEG	1^{st}	LH, CR, CL	
Mn	L	630-670	MEG	1^{st}	LH, CR, CL	
Fe	L	690-730	MEG	1^{st}	LH, CR, CL	SF ₆ for calibration
Со	L	770-810	MEG	1^{st}	LH, CR, CL	
Ni	L	845-885	MEG	1^{st}	LH, CR, CL	
Cu	L	920-950	MEG	1^{st}	LH, CR, CL	
Zn	L	1010-1080	MEG	3 rd	LH	
Ga	L	1105-1165	MEG	3 rd	LH	
Ge						
As	L	1320-1390	MEG	3 rd	LH	
Se	L	1420-1520	MEG	3 rd	LH	
Br						
Kr						
Cd	Μ	390-435	MEG	1^{st}	LH, CR, CL	
Au	Μ	2150-2400	MEG	5^{th}	LH	
Ru	Μ	445-475	MEG	1^{st}	LH, CR, CL	
Ag	Μ	340-390	MEG	1^{st}	LH, CR, CL	
Sn	М	480-510	MEG	1^{st}	LH, CR, CL	

The table is not complete, consult with the Beamline Staff for more information.



Figure B1. (A) Sample Plate in storage case with samples on silicon nitride windows, and labeling information and (B) Sample Holder.



Figure B2. Sample plate coordination system. Note for ALS 5.3.2 Hole #1 and for 11.02 Holes #1,2 and 3 are not accessible.

Appendix B: LOADING/UNLOADING a Sample into the STXM TANK

Appendix B: LOADING/UNLOADING a Sample into the STXM TANK

Recommend that gloves be worn when loading/unloading samples.

Before opening the STXM Tank the following is required:

- 1. Turn OFF the Detector (extents the life of the detector)
- 2. Move the sample holder away from the OSA and Zone Plate by moving ``Coarse Z`` to 5000 (STXM_CONTROL_GUI)
- 3. Close the BSH shutter (on SM_STXM_GUI, protects User from X-rays)
- 4. Open the vent valve on the STXM tank to break the vacuum (may already be at atmosphere)
- 5. Record the sample information, sample position and the empty hole positions (Figs. B1, B2) in your logbook.

Open the STXM Tank by removing the left front plate.

- Remove the sample plate (if applicable) which is held by clips on the Sample Holder (Fig. B1B) by gently lifting it straight up, being careful to avoid hitting the OSA, Zone Plate, detector and mirrors
- 2. Load the sample plate, again being careful to **AVOID** hitting the components in the STXM End Station.
- 3. Look inside the STXM End Station to make sure the sample plate was correctly mounted and the sample holder is about 5000 um back from the OSA. Use a flashlight if necessary.
- **CAUTION**: For **WET CELLS**, pumping may cause the windows to break, thus flushing with Helium may be preferred. In this case do not pump. Ask the Beamline Staff for assistance.
- 4. Make sure that the vacuum pump valve is closed
- 5. Turn on the vacuum pump using the switch on the wall
- 6. Open the vacuum pump valve
- 7. Pump until pressure is less than 200 mTorr (bottom scale on the digital gauge)
- 8. Close the vacuum pump valve
- 9. Turn off the vacuum pump
- 10. Open the Helium valve on the wall (Make sure there is He in the tank, located outside the SM entrance way
- 11. If Helium was not the gas previously used open the release valve to flush the gas line
- 12. Open the gas valve on the STXM tank and allow Helium to fill the tank to 1/6 atmosphere (red line marked on the analog gauge).
- 13. Turn off the Helium valve on the wall
- 14. Turn on the Detector

Appendix C: SM Beamline Polarization Setup

The helical APPLE II type undulator for 10ID-1 Spectromicroscopy beamline is designed to produce plane (Horizontal and vertical), circular, and linear polarized lights with the polarization planes inclined at any angle. It consists of four magnet assemblies that are designed to move longitudinally relative to each other as shown in Figure 7.1. The top and bottom girders can be moved vertically to change the gap and the four assemblies or rails in the longitudinal direction to change the phase.



Figure C1. Schematic of the elliptically polarizing undulator (Attwood 2005).

The following list describes the positions of the magnet assemblies in relation to the operating parameters.

- 1. Circular polarization Two assemblies on the extreme opposite sides are moved in the same direction Parallel mode
- 2. Plane polarization
 - a. Horizontal Four assemblies are at the same position longitudinally
 - b. Vertical Like case 1 for circular polarization, but the longitudinal movement length is equal to the undulator period length
 - c. Inclined Like case 1, but the assemblies are moved in the opposite directions to each other antiparallel mode. i.e., upper left moved forward and the lower right moved backward



Figure D1. Relationship between the position of the zone plate, OSA, sample and detector in the Z direction (e.g., Coarse Z).



Figure D2. Zone plate position as a function of energy. E2 > E1, zone plate moves further away from the OSA.

56

Appendix D: Zone Plate, Order Sorting Aperture (OSA), Sample & Detector Relationship

The spatial relationship between the zone plate, OSA, Sample holder and detector are shown in Figures D1, D2 and D3. OSA is assigned the coordinates (0, 0, 0), left of the OSA is the zone plate and it is negative, while right of the OSA where the sample holder is, is positive. The OSA and Detector do not move during data acquisition, whereas the zone plate and sample holder (i.e., sample) do move during data acquisition.

The zone plate will change position as it is a function of the energy. The Ao value is the expected sample position relative to the OSA. The zone plate uses the Ao in its calculation to determine its relative position to the OSA. Thus, the focal length (f) is calculated based on the distance of the zone plate from the OSA (Z) and the distance from the OSA to the Sample (Ao) (f = Z + Ao). Ao increases with increasing energy. The reason that we change the Ao instead of just moving the zone plate further back from the OSA is due to alignment issues and physical constraints. So if we move from higher energy to lower energy, the calculation of the new focal length is based on the Ao from the higher energy, subsequently the zone plate position is moved closer to the OSA to satisfy f = Z + Ao. Consequently, the zone plate may crash into the OSA when moving from higher energy to lower energy



Figure D3. Schematic of zone plate, OSA and Sample Holder and detector.

Map Network Drive		A Map Network Drive
A	Windows can help you connect to a shared network folder and assign a drive letter to the connection so that you can access the folder using My Computer. Specify the drive letter for the connection and the folder that you want to connect to: Drive: Z: Folder: Image: Connect to a share of the connection and the folder model. Browse Browse Folder: Image: Connect at logon Connect using a different user name. Sign up for online storage or connect to a network server.	What network folder would you like to map? B Specify the drive letter for the connection and the folder that you want to connect to: Drive: Drive: V: Image: Connect to a web site that you can use to store your documents and pictures. Connect to a Web site that you can use to store your documents and pictures.
	< Back Finish Cancel	Finish Cancel

Figure E1. Mapping of Network Drives menus for (A) Windows 7 and B) XP machines.

😫 \\station8107\sm_stxm	<u> </u>
File Edit View Favorites Tools Help	
🛛 🚱 Back 🔻 🕥 👻 🏂 Search 🛛 Folders 🛛 🕼 🔅 🖉 🗶 🎾 🎹 🗸	
Address 😰 \\station8107\sm_stxm	💌 🔁 Go

Figure E2. Windows address window in Explorer

Appendix E: SM Data Access and Data Transfer

The raw data files are archived and saved indefinitely. STXM Data is stored on Drive E. All STXM data is stored by date and the order of collection. For example, the first scan on February 9, 2011 would be stored as 110209000, the next as 110209001, and so forth.

Users can copy their data onto their own laptop or memory stick while at the CLS or can access the data after they have left the CLS via ftp.

Accessing data from the Optical Microscope (VLM Data) or the STXM_Control (STXM Data) computers while at the CLS

Laptops connected to any CLS network port (VLAN110 or 110) can access the internet and transfer data from the STXM_Control or the SM Optical Microscope computer. Wireless is also available, ask the Beamline Staff for the information.

To set-up your laptop requires that you map your network drive and select the appropriate computer and folders.

1. Map the Network Drive.

For Windows 7 computers Click "START" then "COMPUTER" then "Map Network Drive" opening to another screen (Figure E1A).

For XP computers Right Click "START", then select Explore then Right Click "My Network Places" and select "Map Network Drive" opening to another screen (Figure E1B).

2. In the FOLDER option,

For the STXM_Control computer enter <u>\\station8107\sm_stxm\</u>, click Finish For the Optical Microscope computer enter <u>\\WKS-W0000120540\VLM-data</u>\, click Finish

Note that this only needs to be done once on your computer.

Alternatively, can access the data on the STXM_Control computer using the Windows Map_Network_Drive feature by entering '\\station8107\sm_stxm\' in the address window of Explorer (Fig. E2).



Figure E3. Location of USB Extender from the STXM Control Computer.

Memory Stick

STXM DATA: A USB extender on the USB port from the STXM_Control computer makes for easy access (Figure E3). Warning, do not connect the memory stick to the USB port while the STXM is scanning as it may cause the data acquisition to be aborted and/or the computer to crash. The data is stored on Drive E: in the DATA folder

VLM DATA: A USB extender on the USB port from the Optical Microscope computer makes for easy access. The data is stored on Drive C: in the VLM-Data folder.

Accessing data from outside the CLS using an ftp program

In order to access your CLS-SM experimental data from outside the CLS it is necessary that the data be copied to the **ftp data server**, which can only be done from the STXM computer at the CLS. If your data is not in the ftp data server (i.e., not copied yet) and you are no longer at the SM beamline, email the Beamline staff to complete the ftp data transfer.

Data transfer to the CLS ftp data server is done from the STXM_Control computer:

- 1. Double click the "copy STXM data" icon on the desktop. This will copy the data files from STXM_Control to the ftp data server. Note, the VLM data is also accessible provided that it was copied previously from the Optical Microscope using the "copy VLM data" icon on the desktop.
- Use your ftp program to connect to the CLS ftp data server: Port: #21 Username: (ask the Beamline Staff or look in the Manual at the SM beamline) Password: (ask the Beamline Staff or look in the Manual at the SM beamline)

Alternatively type **ftp://ex.lightsource.ca** in internet explorer or in the address line of Windows explorer:

3. Your data is in the \sm-user**STXM-data** folder.

Appendix F: Connecting to the Printer in the STXM hutch

Printer name: prd-cp3525n-e-0

(Windows 7 machines):

- 1. Click "Devices & Printers"
- 2. Click "Add Printer"
- 3. Click "Add Network, wireless or Bluetooth printer"
- 4. Click "The printer that I want wasn't listed" Click "NEXT"
- 5. Select "Select a shared printer by name"
- 6. Enter in box <u>\\vsrv-print-01\prd-cp3525n-e-0</u> (May actually give you list of printers before entering a printer name)
- 7. Click "Next"
- 8. Should indicate you have been successful.
- 9. Click "next"
- 10. Click "finish"

(XP machines)

- 1. Click "Printer and Faxes"
- 2. "File" click "Add Printer"
- 3. Click Next
- 4. Select "A network printer, or a printer attached to another computer" then click "NEXT"
- 5. Select "Connect to this printer (or...)"
- 6. In box type <u>\\vsrv-print-01\prd-cp3525n-e-0</u> (May actually give you list of printers before entering a printer name)
- 7. Then asks you if you want to set as default printer
- 8. Click "Finish"

Appendix G: SM Beamline – Computer's and Software for Users

Users cannot install any software on any SM beamline computer without prior permission from the Beamline Staff.

The program available for image and spectral analysis is aXis2000 (IDL Language), and it is loaded on the STXM Control Computer and the VLM computer.

At the time of data collection (i.e., during your shift), the STXM Control computer can be used for data analysis using aXis2000.

Outside of a User's shift, the SM beamline has one computer (VLM optical microscope) that may be available to Users for data analysis, however, its availability cannot be guaranteed, thus, Users are encouraged to bring their own laptop computers for data analysis.

If you use the SM beamline computers for data analysis, save your work under the folder C:\Users and copy onto a memory stick or transfer through the CLS ftp data server site (see Appendix F) before leaving the CLS as folders older than 1 month will be deleted.

The Beamline staff is familiar with aXis2000 and can assist Users in using it for data analysis. aXis2000 is freeware and is maintained by Adam Hitchcock. The aXis2000 program and the IDL virtual machine package (from ITT Visual Information) can be downloaded from <u>http://unicorn.mcmaster.ca/</u>. Users that do not have an IDL License must use the **IDL Virtual Machine Version**. If you have trouble getting the aXis2000 program to run prior to your arrival at the SM beamline, the Beamline staff will be able to help you install it at the time of your arrival.

1. Optical Microscope (VLM) Computer

(only available when not required for previewing samples) Username: (ask the Beamline Staff or look in the Manual at the SM beamline) Password: (ask the Beamline Staff or look in the Manual at the SM beamline)

aXis2000 – Spectral and image analysis of STXM and PEEM data **Microsoft office suite:** Word, Excel, Powerpoint

Appendix H: Selecting the optimum Cff value

The constant of fixed focus (Cff) value controls the angles for the M1 mirror and the grating of the PGM. Usually the Cff value is 2.15. To change the Cff value consult with Beamline Staff.

One can select different Cff values in order to trade off intensity for increased spectral resolution. This section gives information on how these quantities change as a function of the Cff value for the LEG (250 l/mm) and MEG (500 l/mm).

- Use Cff values between 1.75 to 2.5 for better intensity and to reduce second order light
- Intensity at CirRight polarization is higher than at LinHor polarization.



Figure I1. Peak photon flux in the first order near the carbon, nitrogen and oxygen K-edges in circular polarization measured using the LEG (at STXM exit slit).

Second order contribution near the Carbon K-edge for the LEG and MEG

- 1. Second order contribution increases with an increase in the cff value
- 2. Second order contribution is higher at LinHor polarization than at CirRight polarization
- 3. Use cff values of 2.0 to 2.5 at the Carbon K-edge for both the LEG and MEG.



Figure 8.3. Second order contribution as a function of Cff value near the Carbon K-edge for the LEG (250 l/mm) and MEG (500 l/mm).

10ID-1 CLS SM Beamline Manual

Version: 13 June 2011

Appendix I: References

 Kaznatcheev, K.V., Karunakaran, Ch., Lanke, U.D., Urquhart, S.G., Obst, M., Hitchcock, A.P. 2007. Soft X-ray spectromicroscopy beamline at the CLS: Commissioning results. Nuclear Instruments & Methods in Physics Research A. 582, 96-99.